

# SILICONE BREAST MILK STORAGE SET

## PRODUCT MANUAL

English

**Product Name:** Silicone Breast Milk Storage Set

**Contains:** 2 x 160ml Storage Containers

2 x 250ml Storage Containers

**Materials:** Container Body - Medical Grade Silicone, PPSU

Attachment Ring - PP

Sealing Disk - Medical Grade Silicone

### Features:

- ⦿ Made with medical grade silicone, with great anti-impact performance.
- ⦿ Compatible with Haakaa silicone milk bottle and generation 3 breast pump flange for pumping directly into the bottle. No transfer of breast milk is required.
- ⦿ The twist-tight, leakproof seal makes this container ideal for storing liquids and snacks or taking items on the go.
- ⦿ Air-tight design to retain breast milk's beneficial properties.
- ⦿ Large, easy-to-read millilitre and ounce measurements.
- ⦿ Wide base helps eliminate tip-overs, and the ultra wide-neck design allows for easy cleaning.
- ⦿ Safe in microwave, boiling water, dishwashers and sterilisers.
- ⦿ BPA, PVC and phthalate-free.

### Cleaning & Care:

Wash before and after each use. Although dishwasher safe, we recommend hand washing in warm soapy water and rinsing and drying thoroughly. Do not use any bleach-based cleaners or tablets to sterilise or clean your Haakaa products. To sterilise, use a steam steriliser or boil in water for 2-3 minutes.

### Caution:

- ⦿ Check product condition regularly. Replace if product shows any signs of damage.
- ⦿ Do not store near any sharp objects.
- ⦿ Only soft bristle brushes or soft sponges should be used to clean this product. Hard scourers may scratch the surface.
- ⦿ This product is not a toy. Do not use for any other purpose.
- ⦿ Store product in a cool dry place and avoid direct sunlight.

**Please note:** Volume measurements are an approximation and should be used for general reference only.